

Grillo Gait Trainer

Grillo easily adjusts and is tool-free for user growth, transportation or multi-patient use. Measurement bars allow for quick, easy and accurate adjustment.

Huge growth range: Each frame covers 2 sizes

CONSTANT CENTRE OF GRAVITY with patented height adjustment, allows for postural correction with the user in situ

HUGE GROWTH RANGE means Grillo will accompany the user's growth for a long time within a single size

PRECISE & CONTINUOUS adjustments. Versatile with a range of interchangeable components.

BASIC OR COMPLEX configurations

ANTERIOR **POSTERIOR**

Switch frame configurations if user needs change by reversing the wheels


Direction Walking →



Standard Packages

ANTERIOR PACKAGE CONTAINS

- Pelvic and trunk support
- D Handle
- Pelvic harness
- Front direction locks
- Drag brakes
- Anti-reverse locks
- Foot activated brakes.



POSTERIOR PACKAGE CONTAINS

- Pelvic support
- Pelvic harness
- Trunk support
- Adjustable handle bars
- Direction locks
- Drag brakes
- Anti-reverse locks
- Foot activated brakes



SIZE	POSTERIOR CODE	ANTERIOR CODE	PELVIC SUPPORT HEIGHT	TRUNK SUPPORT HEIGHT	SUPPORT CIRCUMFERENCE	USER HEIGHT RANGE	MAX LOAD
0		5750-0010-000	41-63 cm	52-77 cm	53-79 cm	75-90 cm	35 kg
1-2	5750-0151-12B	5750-0051-12B	49-84 cm	69-100 cm	55-80 cm	90-135 cm	40 kg
3-4	5750-0153-34B	5750-0053-34B	54-100 cm	76-125 cm	67-103 cm	110-160 cm	90 kg
5	5750-0155-05B	5750-0055-05B	63-117 cm	81-143 cm	78-115 cm	up to 190 cm	110 kg