

Accessories

**Leg divider with moving leg cuffs**  
 0 #5750-5400-000  
 1-2 #5750-5401-1-2  
 3-4 #5750-5403-3-4  
 5 #5750-5405-005

**Padded Saddle Seat for Anterior or Posterior**  
 1-2 #5750-0211-1-2  
 3-4 #5750-0213-3-4  
 5 #5750-0215-005



**Arm Gutter Supports**  
 Ant, 1-2 #5750-5501-1-2  
 Ant, 3-4/5 #5750-5503-3-4  
 Post, 1-2 #5750-5511-1-2  
 Post, 3-4/5 #5750-5513-3-4

**Handlebars for Posterior**  
 1-2 #5750-5711-1-2  
 3-4 #5750-5713-3-4  
**Hand grips 1** #5750-5721-1-2  
**Hand grips 3** #5750-5723-3-4



**Height adjust, folding Seat for Posterior**  
 1-2 #5750-0201-1-2  
 3-4 #5750-0203-3-4

**Lumbar Support for Posterior Frame**  
 1-2 #5750-2501-1-2  
 3-4 #5750-2503-3-4



**Weight Bars**  
 1-2 #5750-8431-1-2  
 3-4 #5750-8433-3-4

**Assistant Push Handle**  
 0 #5750-8590-000  
 1-5 #5750-8600-000



**Grab Rail**  
 0 #5750-5700-000

**Parking Brake Pedals**  
 0 #5750-8440-000  
 1-5 #5750-8590-000



Grillo Gait Trainer



**REAR WHEELS FEATURE:**

**Antireverse Device:** Stops wheel from going backward.

**Adjustable Friction:** To reduce speed of movement. 4 levels of adjustment.

**Hand brake compatibility**

**Wheel stop brake:** useful while positioning the child and during transportation.

**FRONT SWIVELLING WHEELS FEATURE:**

**Directional locks:** prevent front wheels from steering, ensure linear stability and allow the user to focus their residual abilities only on linear walking rather than on lateral control.



**QUICK ADJUSTMENTS**

Adjustments are made with levers or handles without need for tools, and can be made with the user in the frame.



**GRADUATED BARS**

On all adjustable supports, to allow reset of the exact position after Grillo has been folded or used.



**SIMPLE & QUICK FOLDING**

Grillo folds without requiring any components to be removed.



**INDEPENDENT TRUNK / PELVIS SUPPORTS**

Ergonomic and fully adjustable in width/ circumference, height, tilt, and forwards/back.



**ERGONOMIC HARNESS**

Has a slight abduction effect without generating pressure. Adjustable, removable and machine washable.



**EASILY STORED AND TRANSPORTED**

Grillo folds in a single operation, reducing size for storage and transportation.