



Atom by Jenx

Technical Specifications



	Measurement	Units	Atom AT01
IMPORTANT INFORMATION: Please be aware, some minimum and maximum measurements are achieved by the use of additional accessory components.			
Age Range	approx.	months - years	6mths - 5yrs
User Weight	max	kg	25
Back Angle	posterior - anterior	degrees	15° - 15°
Back Pad Height from Seat to Top	min - max	mm	270 - 380
Chest Width	min - max	mm	200 - 280
Hip Width	min - max	mm	180 - 270
Seat Depth	min - max	mm	90 - 340
Seat Depth With Long-leg Sitting Board	min - max	mm	315 - 565
Lower Leg Length with Fixed Angle Footplate	min - max	mm	165 - 325
Lower Leg Length with Flip-up Footplate	min - max	mm	115 - 360
Footplate Angle	dorsi - plantar	degrees	5° - 35°
X and Y Base Tilt-In-Space Angle	posterior - anterior	degrees	30° - 10°
Wooden Base Tilt-In-Space Angle High Base and Low Base	posterior	degrees	20°
X Base Footprint	width x length	mm	610 x 620
Y Base Footprint	width x length	mm	650 x 780
Wooden Base Footprint	width x length	mm	500 x 600
X Base Seat Height From Floor	min - max	mm	340 - 510
Y Base Seat Height From Floor	min - max	mm	285 - 750
Wooden Base Seat Height From Floor	high base / low base	mm	550 / 400