

Car Seat Comparisons



Specifications

MILLENIA SICT (AHR)

POSTURE-LUXY

DUALFIX

STARLIGHT-REHA

SPIRIT

IPAI

CARROT 3000

CARROT XL

KIDSLOVER

TIMY

Size 1 Size 2 Size 3 Junior Standard Senior Maxi

DESCRIPTION	Soft approach to early intervention. Postural adaptations behind covers.	One seat, 3 growth configurations: Rear face, forward face and booster modes.	Swivel base makes it easy for parents. Includes posture enhancement kit.	Excellent posture support options and tilt and swivel functions.	Solid shell car seat for larger users. Also suitable for hip spikas.	Fully adjustable special needs car seat with swivel base.	Extension accessories and multi-adjustable features allow longer use.	Larger size Carrot car seat for taller and bigger clients.	Fully adjustable specialised seating with excellent posture support options. Wide range of options and adjustments, allowing complete comfort and postural management.	Correct anatomical seating system which incorporates a pelvic well; often essential to control the pelvis. Offers enhanced comfort and upright sitting.
USER WEIGHT RANGE	Use shoulder height	Use shoulder height	18 kg	9-18 kg, 15-36 kg	59 kg	36 kg	15-36 kg	36-75 kg	22 kg 36 kg 75 kg	15-36 kg 15-36 kg 15-36 kg 15-49 kg
SEAT DEPTH	25 cm	~ 32 cm	20-25 cm	25-30 cm	30 cm or 41 cm	29-42 cm	25-40 cm	45 cm	28-33 cm 34-39 cm 40-46 cm	23-32 cm 32-40 cm 34-42 cm 36-45 cm
SEAT WIDTH	30 cm (Knees)	~ 32 cm	18-25 cm	20-28 cm	38 cm rear, 43 cm front	17-29 cm	31 cm	48 cm	27-31 cm 31-35 cm 36-40 cm	26 cm 31 cm 36 cm 36 cm
SHOULDER HEIGHT	20-40 cm	40 / 46 cm	20-36 cm	32-44 cm	30-50 cm	35-52 cm	32-67 cm	32-67 cm	40-48 cm 48-55 cm 55-64 cm	28-48 cm 35-55 cm 38-58 cm 40-60 cm
APPROX PRODUCT WEIGHT	8.9 kg	~ 7 kg	18 kg	13.5 kg	12 kg	14 kg	8.7 kg	10.5 kg	8.2 kg 8.9 kg 12.3 kg	11 kg 13 kg 15 kg 15 kg
TILT IN SPACE	✓	✓	✓	✓	Tilt bar accessory	✓	Using Wedge	Using Wedge	✓	✓
SWIVEL BASE	✗	✗	✓	✓	✗	✓	Option - requires Isofix	✗	✓	✓
GROWTH ADJUST SEAT	Headrest only	Headrest & Infills	Headrest	✓	✓	✓	✓	✓	✓	✓
FOOTPLATE	✗	✗	✗	✓	✗	✓	✓	✗	✓	✓
TRAY/FOREARM SUPPORT	✗	✗	✗	✓	✗	✓	✓	✓	✓	✓
POSTURAL SUPPORT LEVEL	Mod/Complex	Mild/Mod/Complex	Mild/Mod/Complex	Mild/Mod/Complex	Mild/Mod/Complex	Mild/Mod	Mild/Mod/Complex	Mild/Mod/Complex	Mild/Mod/Complex	Mod/Complex
ISOFIX	Optional	✗	Essential	Optional	✗	Optional	Req. for Swivel	Optional	Optional	✗
REAR FACING CAPACITY	Until child is approximately 2-3 years old	Until child is approximately 6 months old	18 kg	N/A	✗	✗	N/A	N/A	✗	N/A
SUPPORT/PAD OPTIONS	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
TETHER REQUIREMENT?	✓	✓	✗	✗	Optional*	✗	Optional	Optional	✗	Not required*
ANTI-ESCAPE OPTIONS?	✗	✗	✗	✓	✓	✓	✓	✓	✓	✗
RECLINE	✗	✗	✗	✗	✗	✓	✓	✓	✓	Optional Upgrade

*Tether strap may be required for stability or if vehicle design is such that additional restraint is required.

Medifab is as diligent as possible in compiling this information. Medifab does not guarantee the correctness and completeness of the information provided. Equally, Medifab does not guarantee that this information is up to date. For questions concerning our products, their specifications and their utilisation please refer to the instruction manual and/or a product consultant of Medifab. Product specifications are subject to change without notice.



WWW.MEDIFAB.CO.NZ
CALL 0800 543 343
SOLUTIONS@MEDIFAB.COM

SEE PG 51-62 OF YOUR
MEDIFAB CATALOGUE
FOR FULL DETAILS