




Seahorse Plus Toileting & Shower Chair

Shaping better lives.



Feature	Ref	Small Seahorse Plus N72012				Medium Seahorse Plus N72024				Large Seahorse Plus N72036									
		Small Thick Liner - Size 1		Small Thin Liner - Size 2		No Liner		Medium Thick Liner - Size 3		Medium Thin Liner - Size 4		Large Thick Liner - Size 5		Large Thin Liner - Size 6		No Liner			
Seat Depth	SD	327mm	12¾"	355mm	14"	374mm	14¾"	425mm	16¾"	457mm	18"	470mm	18½"	425mm	16¾"	455mm	18"	488mm	19¼"
Rear Width of Seat	SWR	258mm	10"	294mm	11½"	332mm	13"	303mm	12"	339mm	13½"	378mm	15"	411mm	16"	447mm	17½"	477mm	19"
Front Width of Seat	SWF	294mm	11½"	334mm	13¾"	365mm	14¾"	353mm	14"	393mm	15½"	423mm	16¾"	460mm	18"	500mm	19¾"	530mm	20¾"
Height of Back	BH	391mm	15½"	391mm	15½"	414mm	16¾"	490mm	19¼"	490mm	19¼"	514mm	20¾"	549mm	21½"	549mm	21½"	572mm	22½"
Min Headrest height	HH	456mm	18"	456mm	18"	475mm	18½"	562mm	22"	562mm	22"	584mm	23"	620mm	24½"	620mm	24½"	641mm	25"
Max Headrest height	HH	631mm	25"	631mm	25"	650mm	25½"	737mm	29"	737mm	29"	759mm	30"	795mm	31½"	795mm	31½"	816mm	32"
Min height of Seat	SH	540mm	21½"	540mm	21½"	521mm	20½"	521mm	20½"	521mm	20½"	501mm	19½"	521mm	20½"	521mm	20½"	501mm	19½"
Max height of Seat	SH	648mm	25½"	648mm	25½"	629mm	24¾"	625mm	24½"	625mm	24½"	609mm	24"	629mm	24¾"	629mm	24¾"	609mm	24"
Min height of Footrests	FH	275mm	10¾"	275mm	10¾"	255mm	10"	325mm	12¾"	325mm	12¾"	305mm	12"	325mm	12¾"	325mm	12¾"	305mm	12"
Max height of Footrests	FH	375mm	14¾"	375mm	14¾"	355mm	14"	425mm	16¾"	425mm	16¾"	405mm	16"	450mm	17¾"	450mm	17¾"	430mm	17"
Total width (incl anterior support)	OW	532mm	21"	532mm	21"	532mm	21"	586mm	23"	586mm	23"	586mm	23"	692mm	27"	692mm	27"	692mm	27"
Total depth (upright)	OD	807mm	31¾"	807mm	31¾"	807mm	31¾"	860mm	33¾"	860mm	33¾"	860mm	33¾"	877mm	34½"	877mm	34½"	877mm	34½"
Max turning circle (tilted)	TC	1209mm	47½"	1209mm	47½"	1209mm	47½"	1413mm	55¾"	1413mm	55¾"	1413mm	55¾"	1488mm	58½"	1488mm	58½"	1488mm	58½"
Min height over toilet	TH	461mm	18"	461mm	18"	461mm	18"	442mm	17½"	442mm	17½"	442mm	17½"	442mm	17½"	442mm	17½"	442mm	17½"
Max height over toilet	TH	569mm	22½"	569mm	22½"	569mm	22½"	551mm	21¾"	551mm	21¾"	551mm	21¾"	550mm	21¾"	550mm	21¾"	550mm	21¾"
Min width over toilet	TW	426mm	16¾"	426mm	16¾"	426mm	16¾"	478mm	18¾"	478mm	18¾"	478mm	18¾"	592mm	23¾"	592mm	23¾"	592mm	23¾"
Height of Anterior Support to Seat	AH	179mm	7"	179mm	7"	199mm	8"	214mm	8½"	214mm	8½"	234mm	9"	214	8½"	214mm	8½"	234mm	9"
Min Seat angle	SA	0°		0°		0°		0°		0°		0°		0°		0°		0°	
Max Seat angle	SA	40°		40°		40°		40°		40°		40°		40°		40°		40°	
Max User Weight		63.5kg	10st	63.5kg	10st	63.5kg	10st	88.9kg	14st	88.9kg	14st	88.9kg	14st	102kg	16st	102kg	16st	102kg	16st
Product Weight (incl. liners)		21kg	46.3lbs	20.1kg	44.3lbs	18.7kg	41.2lbs	26kg	57.3lbs	24.5kg	54lbs	22.3kg	49.2lbs	29.2kg	64.4lbs	27.7kg	61.1lbs	24.8kg	54.7lbs