

# Walking Frame Comparisons



## Specifications

### GRILLO

### FLUX WALKER

	Size 0	Size 1-2	Size 3-4	Size 5	Size 1	Size 2	Size 3	Size 4
DESCRIPTION	Versatile aid that can be either posterior or anterior ** Interchangeable supports make it suitable for wide range of users.				Lightweight yet sturdy walking aid with unique design that promotes natural walking patterns.			
AGE OF USER*	1-6 years	3-8 years	8-18 years	Adults	1-6 years	3-8 years	6-16 years	Teens/Adults
MAX USER WEIGHT (kg)	35	40	90	110	60	75	75	100
TOTAL BODY HEIGHT (cm)	75-90	90-135	110-160	190	65-115	95-130	110-170	Adult
FOOT - AXILA (cm)	52-77	69-100	76-125	81-143	N/A			
FOOT - HIP (cm)	41-63	49-84	54-100	63-117	N/A			
WEIGHT (Base Package, kg)	14	17	22	31	6.5	6.8	7.1	9.1
ANTI-REVERSE MECHANISM	✓ (Standard feature)				✓ (Standard feature)			
CHEST SUPPORT	Tray	✓	✓	✓	✗			
HIP POSITIONER	✓				✓ (Option)			
FOREARM SUPPORTS	✗	✓	✓	✓	✓ (Option)			
LEG DIVIDER/SEPARATOR	✓				✗			
TOTAL WIDTH (cm)	57 x 72	65 x 75	76 x 91	80 x 98	64 x 65	69 x 81	74 x 88	78 x 91
FOOT TO HAND GRIP HEIGHT (cm)	N/A	52-74	57-90	60-100	48-61(std.) 43-68(uni.)	57-81(std.) 51-88(uni.)	69-87(std.) 62-94(uni.)	82-102(std.) 75-102(uni.)
SADDLE SEAT	✗	✓	✓	✓	✗			
SLING SEAT	✓				✓			
ATTENDANT HANDLE	✓				✗			
FRICTION/SPEED BRAKE	✓ (Standard feature)				✓			
FOLDING SEAT OPTION	✗	✓ Posterior Only	✓ Posterior Only	✓ Posterior Only	✓ (Folding feature can be inhibited at low setting)			

\* Indicator only, for accurate product selection client measurements are the best guide.

\*\* Size 0 only available in anterior configuration

Medifab is as diligent as possible in compiling this information. Medifab does not guarantee the correctness and completeness of the information provided. Equally, Medifab does not guarantee that this information is up to date. For questions concerning our products, their specifications and their utilisation please refer to the instruction manual and/or a product consultant of Medifab. Product specifications are subject to change without notice.