



Shaping better lives.

IRIS Wheelchair Quick-Script

1. PROVIDE USER DETAILS

ASSESSOR NAME		CLIENT NAME			
PHONE		FUNDING	MOH		ACC
EMAIL		TRIAL APPROVAL #			
DELIVERY ADDRESS					

2. WHEELCHAIR TYPE

- ZIPPIE IRIS RIGID
 ZIPPIE IRIS FOLDING
 QUICKIE IRIS

3. SELECT SIZE

W/C WIDTH	MEDIUM PACKAGE	LONG PACKAGE
10"	<input type="checkbox"/>	<input type="checkbox"/>
12"	<input type="checkbox"/>	<input type="checkbox"/>
14"	<input type="checkbox"/>	<input type="checkbox"/>
16"	<input type="checkbox"/>	<input type="checkbox"/>
18"	<input type="checkbox"/>	<input type="checkbox"/>
20"	<input type="checkbox"/>	<input type="checkbox"/>



MEDIUM PACKAGE INCLUDES:


- True-fit growth system
- Tie downs
- Width adjusting system
- Medium seat frame
- 6" x 1.5" semi pneumatic front castors
- 16" Rear mag wheels; airless insert
- 70° legrest hangers with angle adjustable footplates
- Single post height adjustable arm supports with desk length pads
- Height adjustable push handle
- Growable seat pan
- Anti-tips
- Black Opal frame colour
- Foot operated tilt
- Foot operated brakes
- Medium seat depth range: 16-20"

LONG PACKAGE INCLUDES:

- True-fit growth system
- Tie downs
- Width adjusting system
- Long seat frame
- 6" x 1.5" semi pneumatic front castors
- 16" Rear mag wheels; airless insert
- 70° legrest hangers with angle adjustable footplates
- Single post height adjustable arm supports with desk length pads
- Height adjustable push handle
- Growable seat pan
- Anti-tips
- Black Opal frame colour
- Foot operated tilt
- Foot operated brakes
- Long seat depth range: 18-22"

4. SET SPECIFICATIONS

SEAT WIDTH (R)		SEAT TO FLOOR	
SEAT DEPTH (S)		HANGAR ANGLE	
SEAT TO FOOTPLATE		BACK HEIGHT (T)	



5. CONFIGURE OPTIONS

ARM OPTIONS	<input type="checkbox"/> Full length arm pads	<input type="checkbox"/> JAY padded arm trough	<input type="checkbox"/> Tray
BACK OPTIONS	<input type="checkbox"/> Reclining backrest	<input type="checkbox"/> Extra tall back canes	<input type="checkbox"/> Black top <input type="checkbox"/> Clear top
HARNESS OPTIONS	<input type="checkbox"/> Centrepoint harness	<input type="checkbox"/> Vest harness	
HIP BELT OPTIONS	<input type="checkbox"/> 2 Point hip belt	<input type="checkbox"/> 4 Point hip belt	
LEG OPTIONS	<input type="checkbox"/> Elevating legrests	<input type="checkbox"/> 80° legrest hangers	<input type="checkbox"/> Calf strap
FOOT OPTIONS	<input type="checkbox"/> Contracture footplates	<input type="checkbox"/> JAY Footbox	

6. NOTES