

Quickie 2 Wheelchair Quick-Script



1. PROVIDE USER DETAILS

ASSESSOR NAME		CLIENT NAME			
PHONE		FUNDING	MOH		ACC
EMAIL		TRIAL APPROVAL #			
DELIVERY ADDRESS					

2. SELECT SIZE - QUICKIE 2

Standard stocked sizes

		W/C DEPTH			
		16"	17"	18"	20"
W/C WIDTH	16"	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>
	17"		<input type="checkbox"/>		
	18"	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>
	20"			<input type="checkbox"/>	<input type="checkbox"/>

Standard package in other size ____ x ____ (Indent order)



PACKAGE INCLUDES:

- DAB depth & angle adjustable backrest
- 24" lite spoke rear wheels
- Anodised push rims
- Marathon plus tyres
- 5" x 1.5" semi-pneumatic castors
- Single post height adjustable armrests with full length pads
- 80 degree footrest hangers
- Angle adjustable footplates
- 16" high 3DX tension adjustable back upholstery
- Integrated push handles
- Anti-tips
- Spex 2 point belt
- Colour is Black Opal

3. SET SPECIFICATIONS

SEAT WIDTH (R)		CAMBER		(3°)
SEAT DEPTH (S)		COG		(2.5")
FRONT SEAT TO FLOOR	(18")	SEAT TO FOOTPLATE		
REAR SEAT TO FLOOR	(17")	BACK ANGLE		(90°)
	(Default)			(Default)

4. STANDARD CONFIGURATION OPTIONS

<input type="checkbox"/> 70 degree extn mount hangers	<input type="checkbox"/> Desk length armrests	<input type="checkbox"/> ABS Side guards instead of Single Post armrests	<input type="checkbox"/> Neoprene impact guards
<input type="checkbox"/> 4"x1.5" castors	<input type="checkbox"/> 6"x1.5" castors	<input type="checkbox"/> Spoke Guards 24"	

5. NOTES