

# REAL 9000 Modular Work Chair

Maintain your independence! REAL 9000 helps users to be more independent by enabling them to perform tasks without assistance either at home or in a working environment. An easily adjustable seat height (electrically or manually) ensures an ergonomically correct position at the work surface.

*Supreme functionality, quality and ergonomics*

#3760-0010-000

**EXTRA SUPPORT** at the sides and between the shoulder blades

Readily accessible, **ERGONOMICALLY DESIGNED CONTROLS** which can be positioned on the right or the left by the seat or under an armrest

Switch between **ACTIVE & PASSIVE** position with seat angle adjustment

**EASY TO CONTROL**, width-adjustable brake lever bracket, fitted on right and left side

**COMPACT MOBILE BASE** offers plenty of space for your legs and feet



Move around home or work with castors that turn smoothly in all directions.

Angle adjust to allow position change over the course of the day to relieve spinal pressure.

Increase reaching range while seating with seat height adjustment of 28 cm (electric model)

## Designed for people with:

- Post injuries; perform tasks with ease again while recovering from bodily injuries
- Specialist seating requirements and physical disabilities
- Elderly care; Mobility around home and kitchen

## Models:

### REAL 9000 Plus:

Select manual or electric height adjustment, and standard or Coxit seat (divided with individually adjustable flaps to allow the user to sit with an open hip angle).



### REAL 9000 Plus Child:

Paediatric version of the REAL 9000 to accommodate a child's needs. The ergomedic plus seat system and a large range of functional accessories are available as optional extras.



### REAL 9200 Twin Electric Plus:

Modular, bariatric chair with additional reinforcement, ideal as a mobility aid for slightly heavier users. Able to withstand weights up to 275 kg.



### REAL 6100 Plus

Electrically powered wheelchair for indoor use. Seat easily moves on small surfaces and thresholds in the home, at work, or in school.



REAL CHAIR	MOBILE BASE	SEAT CUSHION SIZES	BACKREST SIZES	SEAT HEIGHT	SEAT ANGLE	USER WEIGHT
SPECS IN CM	48 x 53	Adult & Coxit: 44 x 48, 48 x 48, 44 x 43. Adult: 40 x 43, 40 x 39, 40 x 48	Low 33 x 26 Low 38 x 26 High 33 x 36 High 40 x 34	All intervals can be adjusted +6 /-2 cm Std manual: 45-65. Std electric: 47-75. Low Manual: 40-54. Low Electric: 40-60	+15/ -8°. Crank adjust; Gas spring, seat controls; Gas spring, armrest controls; Electric. Static option.	(Specify when ordering) <55 kg, 55-90 kg, 90-130 kg manual, 90-150 kg electric, >130 kg (custom)

REAL 9000