

Spaghetti Shower Chair

Unique backrest design has been proven to work when others fail! Flexi-Tube upholstery assumes the shape of the user and provides exceptional support during showering, even including individuals with scoliosis and other asymmetries.

Comfort with simplicity for many body shapes



High Wear - resistant PVC Tubing

Adjustable headrest (included in standard package)

Light and compact for carer-friendly showering

15° seat dump angle assists with ease of hygiene during showering

Durable 12cm wheels with rear brakes




Removable armrests for ease of transfers

Elevating, height adjustable legrest to accommodate reclined postures

Triangular shaped backrest with recline holds the user comfortably and securely without need for laterals.

FLEXI-TUBE UPHOLSTERY

- Upholstery warms and softens when being used, therefore moulding to the body's shape.
- This means longer tolerance is maintained during showering.
- Provides comfort for mild to severe disabilities and body shapes.
- Save spending large amounts on customising to individual contours.



ACCESSORIES

Abduction Pommel
6071-1500-210




Padded Commode Seat*
6071-0419-812



Tubular Padded Arm Support
6071-5550-049



Padded Lap Belt
6071-6203-017



Comfort Seat Liner
6071-0450-209



Comfort Back Support Liner
6071-2450-209



Code	Age Range	Max User Weight	Seat Depth	Seat Width	Seat To Floor Height	Footplate Depth Adj.	Recline Angle	Seat Dump Angle	Chair Weight
6071-0000-210	Adolescent to Adult	110 kg	40 - 46 cm	46 cm	Front: 57 cm Rear: 45 cm	36 - 46cm	90°-115°	15°	18 kg