

Starlight Car Seat

Fully adjustable special needs seat, offering excellent posture support options and an inbuilt tilt function. Starlight is suitable for young children from 9kg, and will grow with the child right up to 36kg.

Excellent Postural Support




- PACKAGE INCLUDES:**
- Backrest infill padding
 - Soft shoulder padding
 - Slim lateral trunk supports
 - 10° tilt
 - Base board accessory
 - Swivel base system (left turning)
 - Seat stabilisation clamp bow
 - Grey upholstery



- SLIM SWIVEL BASE**
Allows more space and headroom inside the car when turning out.
- MANY OPTIONS**
For comfort, functional and postural requirements.
- BACKREST RECLINE**
With inbuilt one position tilt function at 10°.

Accessories

- ISOfix adapter**
For stability in vehicle
#4140-6910-100
- Seat Protection Mat**
Protects vehicle upholstery
#4140-0190-000



- Adjustable Footrest**
#4140-5240-000




- Padded Foot Roll**
#4140-5340-000



- Forearm Support Tray**
#4140-5530-000



- Abduction Pommel**
#4140-1510-000



- Longer Harness Straps**
Additional 15 cm
#4140-6570-000



- Extra 15° Tilt system**
#4140-0241-015



CODE	USER WEIGHT	SEAT WIDTH	SEAT DEPTH	SHOULDER HEIGHT	CAR SEAT WEIGHT
4140-0011-090	9-36 kg	20-28 cm	25-30 cm	32-44 cm	13.5 kg