

# Vehicle Seating Comparison Chart



SPECIFICATIONS	HERO		STARLIGHT	CARROT 3000	CARROT XL	KIDSFLEX		
	With Standard Cushion	W/out Standard Cushion				Size 1	Size 2	Size 3
DESCRIPTION	Swivelling special needs car seat for forward or rear facing.		Excellent posture support option with tilt and swivel functions.	Extension accessories and multi-adjustable features allow longer use.	Larger version of Carrot 3000 for older or heavier children and small adults.	Fully adjustable specialised seating with excellent posture support options. Available with a wide range of options and adjustments, allowing complete comfort and postural management for the child.		
USER WEIGHT RANGE (kg)	18 (Max.)		9-36	15-36	36-75	22	32	75
SEAT DEPTH (cm)	25	28	30	25-40	38 or 42	28-33	34-39	40-46
SEAT WIDTH (cm)	21	25	30	31	48	27-31	31-35	36-40
SHOULDER HEIGHT (cm)	21-34	27-35	27-37	32-67	50 - 65	40-48	48-55	55-6
APPROX PRODUCT WEIGHT (kg)	16		13.5	8.7	10.5	8.2	8.9	12.3
TILT IN SPACE	✓		✓	Using wedge	Using wedge	✓		
SWIVEL BASE	✓		✓	Option - requires ISOFIX latch attachment	✗	✓		
GROWTH ADJUST SEAT	✓		✓	✓	✓	✓		
FOOTPLATE	✗		✓	✓	✗	✓		
TRAY/FOREARM SUPPORT	✓		✓	✓	✓	✓		
POSTURAL SUPPORT LEVEL	Mild/Mod/Complex		Mild/Mod/Complex	Mild/Mod/Complex	Mild/Mod/Complex	Mild/Mod/Complex		
IS ISOFIX REQ. IN VEHICLE?	Essential		Optional	Only req. for Swivel Accessory	Optional	Optional		
REAR FACING CAPACITY	✓		N/A	N/A	N/A	✗		
SUPPORT/PAD OPTIONS	✓		✓	✓	✓	✓		
TETHER REQUIREMENT	✗		Seatbelt must be used behind client until 18 kg		Optional	Optional		✗
ANTI-ESCAPE OPTIONS	✗		✓ Custom		✓	✓		✓ Custom
RECLINE	✗		✗		✓	✓		✓

\* Indicator only, for accurate product selection client measurements are the best guide.  
\*\* With flat cushion

Medifab is as diligent as possible in compiling this information. Medifab does not guarantee the correctness and completeness of the information provided. Equally, Medifab does not guarantee that this information is up to date. For questions concerning our products, their specifications and their utilisation please refer to the instruction manual and/or a product consultant of Medifab. Product specifications are subject to change without notice.