

# Indoor Seating Comparison Chart



SPECIFICATIONS	JENX CORNER SEAT			COBRA BASE			JENX ATOM	JENX JUNIOR+ SEAT			JENX MULTISEAT		GRAVITY CHAIR				KRABAT JOCKEY LITE			KRABAT JOCKEY PLUS		
	Size 1	Size 2	Size 3	Size 0	Size 1	Size 2		Size 1	Size 2	Size 3	Size 4	Size 5	Sml	Med	Lrg	XLrg	Size 1	Size 2	Size 3	Size 1	Size 2	
DESCRIPTION	Popular as the child's first seating system at floor level.			High & Low level seating for the Bingo or Spex Wonderseat			A modular seating system for users 6 months - 5 years	Comfortable, contoured and stable modular seating system, with outstanding levels of adjustment and a wide variety of support options.			Larger version of the Jenx Junior+ Seat, designed to suit older children and adolescents. Indent order only (non-stocked).		Supportive alternative seating and foam positioning system to augment postural control for a wide range of individuals. Also used as a 'hugging' chair for individuals who benefit from this type of sensory input.				The Krabat Jockey Lite is the newest addition to the series of active children's chairs from Krabat. The compact and modern design of the undercarriage makes it easy to use even where space is limited. The chair is height adjustable and the seat can be tilted.			The Krabat Jockey Plus is an active chair with a saddle seat that offers more support and a higher level of comfort than the other Jockey chairs. The Plus backrest supports the upper body with adjustable side supports, head support and the use of different harnesses.		
APPROX AGE OF USER*	6 mths-2 yrs	1-6 yrs	4-10 yrs	4 mths-4 yrs	1-5 yrs	3-15 yrs	6 mths - 5 yrs	18 mths - 5 yrs	3-8 yrs	5-12 yrs	8-14 yrs	10-18+ yrs	Refer to Measurements				2-7 yrs	7-14 yrs	14 yrs+	2-6	6-14	
SEAT DEPTH (cm)	35	37	46	12-27	16-34	22-45	9-34	21-30.5	24-32	27.5-39	28.5-45	35.5-52	26	32	40	45	28	34	42	26	34	
SEAT WIDTH (cm)	33	39	49	13.25	18-30	23-35	18-27		15-37.5		20-44		27 front 20 cm back	30.5 cm front 23 cm back	37 cm front 27 cm back	40 cm front 30 cm back						
SHOULDER HEIGHT (cm)	31-42	44-62	50-68	25.5-39.5	28-46	34-63	27-38	31-41	33-43	39-52	46.5-59	52.5-66	60cm total back	69cm total back	74cm total back	80cm total back						
FOOTPLATE HEIGHT RANGE (cm)	N/A			16-34	16-37	16-40	16.5-32.5 (Fixed Angle Footplate) 11.5-36 (Flip Up Footplate)	19.5-27	25-36.5	25-40.5	25-47		24 cm to floor	26 cm to floor	28 cm to floor	30 cm to floor	14-23	14-23	14-30			
MAX LOAD (kg)	25 kg	35 kg	55 kg	25	50		25		50		100		50	70	100	130	80	80	100	Manual: 65 Electric: 50	Manual: 65 Electric: 50	
APPROX PRODUCT WEIGHT (kg)	Approx. 16			Seat: 7.5 Buggy: 17.5 (total)	11.3	11.5	13	~20kg (variable)			Variable		4	5.2	7.4	9.6	10	10	11.6	Manual: 12.2 Electric: 13.9	Manual: 16 Electric: 17.7	
TILT IN SPACE	✗			45°			✓	-10° - +30°			-10° - +30°		Use wedge underneath				-15° - 0°	-15° - 0°	-15° - 0°	-20° - +20°	-20° - +20°	
RECLINE	✗			80°-180° in 10° steps			90°-160°	-15° - +15°			Adjustable Angle		✗									
GROWTH ADJUST SEAT	✓			✓			✓	✓			✓		✗									
SEAT HEIGHT (cm)	N/A			27-61			28.5-75	29-76	33-76	36.5-76	33-57		N/A				45-57	45-57	57-75	39-79	40-80	
HI-LOW ACTIVATION	N/A			Hand Operated Gas Assist			Hydraulic Foot Pump	Hydraulic Foot Pump			Hydraulic Foot Pump		N/A									
TRAY	✓			✓			✓	✓			✓		✓									
POSTURAL SUPPORT LEVEL	Mild/Moderate			Moderate/Complex			Mild/Moderate/Complex	Mild/Moderate/Complex			Mild/Moderate/Complex		Mild/Moderate/Complex									

\* Indicator only, for accurate product selection client measurements are the best guide.

Medifab is as diligent as possible in compiling this information. Medifab does not guarantee the correctness and completeness of the information provided. Equally, Medifab does not guarantee that this information is up to date. For questions concerning our products, their specifications and their utilisation please refer to the instruction manual and/or a product consultant of Medifab. Product specifications are subject to change without notice.