

# Q-Series Powered Wheelchairs



BASE	Q700			Q500M	Q500M Mini	Q400M	Q300M Mini
	Q700R	Q700M	Q700F				
SEDEO UP		✓	✓				
SEDEO ERGO	✓	✓ C-ME	✓				
SEDEO PRO	✓	✓ C-ME	✓	✓ C-ME	✓		
SEDEO LITE	✓	✓ C-ME	✓	✓ C-ME		✓	✓
CAPTAINS SEAT		✓		✓		✓	✓
DRIVE TYPE	Rear	Mid	Front	Mid	Mid	Mid	Mid
TURNING RADIUS	950 mm	560 - 630 mm	735 mm	560 - 630 mm	550 mm	560 - 630 mm	550 mm
DRIVE WHEELS	14"	14"	14"	13"	12.5"	13"	12.5"
CASTOR WHEELS	9" & 10"	6"	9" & 10"	7"	7"	7"	7"
MAX REVERSING WIDTH	1180 - 1370 mm	1100 mm	1240 mm	1100 mm	950 mm	1100 mm	950 mm
MAX CURB CLIMBING	100 mm w/KC, 80 mm w/o KC	80 mm w/o KC	100 mm w/o KC	50 mm	40 mm	50 mm	40 mm
MIN TOTAL WIDTH	640 mm	622 mm	640 mm	620 mm	520 mm	620 mm	520 mm
TOTAL LENGTH (MAX)	1180 - 1370 mm	1100 - 1190 mm	1340 mm	1280 mm	1020 mm	1280 mm	1120 mm
WEIGHT CAPACITY	160 kg	160 (120 UP) kg	160 (120 UP) kg	136 kg	136 kg	136 kg	136 kg
MAX BATTERY RANGE	40 km	40 km	40 km	40 km	30 km	40 km	30 km
GROUND CLEARANCE	65 mm	70 mm	60 mm	70 mm	80 mm	70 mm	80 mm
APPROX WEIGHT	175 kg 198 kg (UP only)	175 kg 198 kg (UP only)	175 kg	165 kg	103 kg	165 kg	103 kg

AVAILABLE WITH	SEDEO UP	SEDEO ERGO	SEDEO PRO	SEDEO LITE
	Q700 M F	Q700 R M F	Q700 R M F Q500 M Q500 M Mini	Q700 R M F Q500 M Mini Q400 M Q300 M Mini
SEAT FRAME SPECS				
LOWEST SEAT TO FLOOR HEIGHT	450 mm	450 mm	415 mm	410 mm
SEAT WIDTH	400 - 500 mm	400 - 560 mm	410 - 560 mm	385 - 510 mm
SEAT DEPTH	400 - 540 mm	400 - 560 mm	400 - 600 mm	410 - 560 mm
BACKREST HEIGHT	560 - 710 mm	560 - 710 mm	440 - 620 mm	460 - 510 mm
BACK REST ANGLE ADJ	NA	NA	90 - 120 degrees	86 - 106 degrees
ARMREST HEIGHT FLIP BACK	260 - 350 mm	260 - 350 mm	260 - 350 mm	260 - 350 mm
ARMREST HEIGHT SINGLE POST	NA	NA	260 - 360 mm	260 - 360 mm
WEIGHT CAPACITY	120 kg	160 kg	136 kg (Q700 -160 kg)	136 kg
POWER SEAT FUNCTIONS				
ELECTRONICS	R-Net only	R-Net only	R-Net only	VR2 & R-Net
ANTERIOR TILT	25 degrees	NA	NA	NA
TILT	22 degrees	50 degrees	30, 50 degrees	30, 50 degrees
MANUAL RECLINE	NA	NA	90 - 130 degrees	85 - 120 degrees
POWER RECLINE	NA	NA	90 - 130 degrees	85 - 120 degrees
POWERED RECLINE - SHEAR REDUCTION	175 degrees	175 degrees	90 - 150 degrees	NA
LIFT	200 mm	300 mm	300 mm	300 mm

**C-ME** No more waiting around or driving at a reduced speed when elevated. C-Me enables you to elevate up to 12" while traveling up to 5 kph! You can lift or descend at the touch of a button with no need to stop while transitioning. Just rise-up and go at the same time!

## DESIGN AND APPEARANCE

### DISTINGUISHED STYLE

With the Q-Series, we've taken style to the next level. With a matt-black chassis finish, hidden wiring and components combined with sharp-tailored upholstery, it's a sleek and a contemporary design with a variety of colours to choose from.



\*Q700/Q500: ● ● ● ● ● ● ● ●

\*Q500 Mini: ● ● ● ●

\*Q400: ● ● ● ● ●

\*Q300: ● ● ● ●

### FANCY A CHANGE?

A choice of contemporary colour inlays to perfectly complement its striking matt-black bodywork.

Simply replace the colour inlays. Styling your powered wheelchair has never been so simple.



\* The number of colours available may differ depending on the seat and base combination.

## CONTROLS AND CONNECTIVITY

### VR2 & R-NET

#### ADVANCED JOYSTICKS

The QUICKIE Q-Series wheelchairs come with the reliable VR2 control or the superior R-Net advanced joystick. A large, easy-to-read HD display and simple to control paddle switches allow for intuitive operation. (Depends on model).



## ALTERNATIVE CONTROLS

### SWITCH-IT

#### INNOVATIVE CONTROLS

Conveniently access your 5 programmable positions with the Ctrl+5 – a perfect combination with Memory Seating (Memory Seating with SEDEO ERGO and SEDEO ERGO UP only). It allows you to go straight to the function you need at the touch of a button.



### LINK-IT

#### MOUNTING SYSTEM

Need to use additional switches or special controls? No problem, easily connect them using the universal Link-It mounting system.

