

# Ella Stander

A multiplicity of standing options and huge growth range maximises your investment in this versatile aid for individuals with moderate to severe disabilities.

## Versatility for Prolonged Standing



**Highly durable frame:**  
easy to assemble and adjust

**Easy care upholstery**  
can be removed for machine washing

**Stepless support**  
adjustments ensure accurate and comfortable positioning

**Individual foot and knee adjustments**

Can be set up for **supine or prone standing:** change over set-up with only 3 steps

## Knee and ankle flexion adjustment for:

- Angles in all directions (including abduction)
- Height
- Width
- Depth



Prone set-up



Supine set-up (shown without tray table)

Good hoist access and swing-away side supports to facilitate easy transfers for supine standing

## ACCESSORIES

Standard package includes: headrest, tray, body supports and knee cuffs and offset brackets.

**Shoulder Protractors**  
5220-8460-000



**Harness Fitting Kit**  
5220-8460-000



**Swing Away Laterals**  
5220-8460-000



**Tray with Bowl**  
5220-8460-000



## ELECTRIC POWER ASSIST

Assisted 90° electric tilt function makes ease of adjustment, reduces manual handling and improves client experience.



## EASY TO ORDER

Comes standard as fully complete package

Size	Code	Max User Weight	Approx User Height	Footplate to Axilla	Chest / Hip Width	Foot to Petella	Foot to Pelvis	Abduction (Max Angle)	Kneepads - Footrest
1	5221-0020-000	30 kg	60-95 cm	43-68 cm	12-20 cm	14-45 cm	30-54cm	50°: 60 cm User Height 42°: 95 cm User Height	0-25°
2	5222-0020-000	40 kg	90-130 cm	57-97.5 cm	17-25 cm	18-55 cm	47.5-78.5cm	44°	0-25°
3	5223-0020-000	55 kg	120-155 cm	80-128 cm	19-30 cm	20-65 cm	50-104 cm	40°	0-25°
4	5224-0020-000	80 kg	130-180 cm	90-150 cm	25-41 cm	25-85 cm	62.5-122.5cm	40°	0-25°