

Hero

—
SAFETY FEATURES

Swivelling special
needs carseat
for forward or
rear-facing



Hero

Hernik Hero car seat is adapted from the Recaro range and is specifically modified to accommodate for special needs children.

- The Hero is designed to be used in the back seat of a vehicle and **does not require use of the vehicle seat belt; the Isofix anchorage system providing a secure attachment to the vehicle.**
- Set the switch (located on the underside of the car seat) to **rear or forward setting. Rear setting must be used for a child of 0-15 months**, although in NZ recommendation is until 2 years old.
- Install the child seat in the car in such a way that it is **not trapped by the front seat or the vehicle doors.**
- The **support leg must be extended to the appropriate length** to suit the vehicle.
- The **belts must not be twisted or trapped** and must be completely **tightened.**
- **Before each journey, check that the indicators on the front of the seat base are green.** These indicators show that the following is correct: Isofix secured, support leg in place, and swivel locked.
- Ensure child has **not outgrown the weight capacity which is 18kg.** If child is approaching this weight limit, consider what car seat would be next option.
- Periodically **check for loose or missing parts.**



Ensure child has not outgrown the weight capacity limit

Contact the car seat specialists for more information:
0800 543 343
solutions@medifab.com

Distributed by **medifab**[®]